

AS 1

Alternating Pressure System

Prevention and therapy of decubitus

- | for the treatment of decubitus up to stage I
- | Continuously adjustable control unit
- | External filter
- | Micro-perforation ventilation overlay
- | Extra material flaps on both ends



Technical information

Control unit

- | Size: 24 x 12 x 9.5 cm
- | Weight: 1.4 kg
- | Cycle: 10 min.

Mattress

- | Size: 200 x 90 x 6.5 cm
- | Weight: 1.6 kg
- | Material: PVC, anti-flammable
- | Max. patients weight: 90 kg

